

## MEHROTRA BIOTECH RH VIEW BOX



Mehrotra Biotech RH View Box comes with a soft fluorescent, glare free LED that provide excellent uniform slide illumination for easily and accurately monitoring the viewing area.

RH view box is having application in Blood typing, RH determinations, Warming slides for Gram's strain, Reading enzymes antibody screen in micro plates.

### Features:

- Blood Typing.
- RH Determinations.
- Warming slides for Gram' strains.
- Tissue typing.
- Reading enzymes antibody screen in micro plates.

### Technical Specifications:

Construction material	AISI 304 stainless steel with top lid made up of acrylic.
Dimension (D*W*H)	24*24*5inch
Illumination	LED
Switch	On/ Off
Power	220V, 50Hz
Power consumption	30W
Indicator	For Power supply
Weight	3 Kg

ISO 9000: 2015

ISO 13485: 2016

**MEHROTRA BIOTECH PVT. LTD.**

M&M TOWER B-22 23 VIBHUTI KHAND GOMTI NAGAR, Lucknow 226010

Tel + 91 522- 2720663

E-mail: [sales@mehrotrabiotech.com](mailto:sales@mehrotrabiotech.com), [service@mehrotrabiotech.com](mailto:service@mehrotrabiotech.com)

Website: [www.mehrotrabiotech.com](http://www.mehrotrabiotech.com)

