

MEHROTRA GROUP



**HOUSE OF BIOMEDICAL RESEARCH
& SCIENTIFIC EQUIPMENT**

BLOOD PLASMA THAWING MACHINE



MBTM-W



MBTM-D

➤ADUANTAGES

- 1.Thick film instantaneous water and electricity separation heating system,
- 2.Freely adjustable defrost barrier,
- 3.Intelligent operation management system,
- 4.Automatic water filling function,
- 5.USB data recording and storage,
- 6.Compact and flexible, easy to place and use


➤APPLICATION

MBTM series is mainly used for constant temperature heating and thawing of plasma, platelets and blood cells. It is also suitable for constant temperature heating of medical frozen products, blood transfusion bottles(bags), anesthetics and medicinal crystals.

This model is also suitable for constant temperature heating of fiber adhesives, medical frozen products, and operating room anesthetics.

Model	MBTM-W	MBTM-D
Thawing Quantity(Optimal/Maximum)	6 bags only	3/6 bags
Temperature Control Range	RT~60℃	
Temperature Control Accuracy	±0.1℃	
Water Storage Capacity	22kg	
Thawing Time	20~30 minutes	
Rated Consumption	1500W	
Power Supply	AC220/110V,50/60Hz	
Dimensions(WxDxH)(mm)	550x460x450	
Shipping Size(WxDxH)(mm)	670x560x580	
N.W./G.W. (kg)	30/35	

MEHROTRA GROUP



HOUSE OF BIOMEDICAL RESEARCH & SCIENTIFIC EQUIPMENT

M&M TOWER B-22 23, VIBHUTI KHAND, GOMTI NAGAR, LUCKNOW- 226010

Tel + 91 522- 2720663, 9838203524,9838203399,9839014323

E-mail: sales@mehrotrabiotech.com, service@mehrotrabiotech.com.

Website: www.mehrotrabiotech.com

OUR BRANCH OFFICES

LUCKNOW

DELHI

MUMBAI

HYDERABAD

BANGALORE

KOLKATA

AHMEDABAD

