

# TD5G Table type low speed centrifuge



## **Features:**

1. Brushless DC motor, maintenance-free.
2. Controlled by Microprocessor with touching panel and digital display, easy to operate.
3. Steel body, stainless steel centrifugal chamber, resist corrosion, easy to clean, safe and stable.
4. Arbitrarily setting the parameter of time and speed etc, easy to use.
5. 0-9 Accel/Decel modes, 40 self-defining programs can be stored and recall.
6. Short spin function.
7. It's mainly used in biochemistry Laboratory and Clinical Laboratory in Hospital and so on.

## **TECHNICAL PARAMETER**

Model	TD5G
Max Speed	5000r/min
Max RCF	4020xg
Max capacity	300ml
Time Range	0~999min/sec
Speed Accuracy	±20r/min
Noise	≤60db
Power Supply	AC 220V 50Hz/110V 60Hz, 8A
Dimension	340×400×270(mm)
Weight	20kg

## **Rotor Options:**

Rotor Type	Rotor No.	Capacity	Max Speed (r/min)	Max RCF (×g)	Dimension of Tube (mm)
Angle Rotor	No. 38	4×50ml	5000	4020	Φ38×124
	No. 38	6×50ml	5000	4020	Φ38×124
	No. 37	12×10ml	5000	4020	Φ16×82
	No. 37	12×15ml	5000	4020	Φ16×115