









WHAT WE DO?

Mehrotra Biotech Pvt Ltd has been established with an aim to contribute to the enhancement of medical-hospital, pharmaceutical, biotechnological, and general laboratory services through the provision of high-quality innovative solutions, ensuring an improved quality of life and well-being for society at large.

Our primary objective is to deliver customized solutions that precisely align with the unique requirements of our customers. To achieve this, our dedicated functions including engineering, R&D, manufacturing, and customer service are focused on each technology within Mehrotra Biotech's area of expertise:

- We offer solutions and equipment for clean air technology, encompassing both horizontal and vertical laminar flow cabinets, biological safety cabinets, and laminar flow systems.
- Our portfolio includes solutions and equipment for advanced sterilization utilizing Autoclaves, Ovens, and incubators.
- We specialize in solutions and equipment employing biological containment technology, such as containment isolators and sterile isolators.
- Additionally, we provide solutions and equipment utilizing freeze drying and vacuum drying technology, ranging from laboratory freeze dryers to GMP freeze dryers.

By offering these comprehensive solutions, we aim to continually improve and advance the various sectors we serve, ensuring the highest standards of quality, safety, and efficiency in medical and laboratory practices.

"Our aim is to provide tailor-made solutions in strict accordance to the needs of customers."

WATER BATH



Mehrotra Biotech make Water bath WB series provides precise temperature control to create a stable temperature environment for a variety of applications. Easy to use and maintain, leakproof stainless steel chamber available with a digital control system. Used in Ideal for warming reagents, bacteriological examinations, microbiological assays, sample thawing, coliform determinations and environmental and food technology applications.

TECHNICAL SPECIFICATIONS

Model	WB05	WB8	WB12	WB20	WB30
Capacity	5L	8L	12L	20L	30L
Temp Range	RT+5°C to 100°C				
Temp Stability	±1.5°C				
Temp Accuracy	±0.5°C				
MOC(Inner)	Stainless Steel AISI304				
MOC(Inner)	MS Powder Coated				
Concentric Ring (Optional)	-	2 Rows, 4 Holes	2 Rows, 6 Holes	2 Rows, 8 Holes	2 Rows, 10 Holes
Ext Dimension (LxWxH) in	15x12x12	15x14x12	17x 14x13	22x15x13	22x10x15
Sensor	PTS100				
Controller	Dual LED Display programmable controller				
Power	500W	750W	1000W	1200W	1500W
Power Supply	220V, 50/60Hz 220V, 50/60Hz				
*Note: Product is customizable as per user requirement.					

AVAILABLE DESIGN OPTIONS





Mehrotra Biotech Pvt. Ltd.

M&M TOWER B-22 23, VIBHUTI KHAND, GOMTI NAGAR, LUCKNOW- 226010

Tel + 91 522- 2720663, 9838203520

E-mail: sales@mehrotrabiotech.com, service@mehrotrabiotech.com.

Website: www.mehrotrabiotech.com

OUR BRANCH OFFICES

