

SLEEK DONOR COUCH SC100

Minimizes vasovagal attack and
motivate repeat voluntary donors

ISO 13485



Optimize blood care with advanced technology

“ To attract volunteer donors and to encourage their continued participation, it is essential that the conditions surrounding blood donation be as pleasant, safe and convenient as possible”
AABB Manual

“ Fainting : Place the donor on his/her back and raise the feet above the level of the donor’s head”



FEATURES

- **Best suitable** for mobile blood donation van and blood banks.
- Ensures **safety** and comfort to the donor and phlebotomist.
- Single **switch** operation in the instance of **Vasovagal Attack**.
- The **soft cushioned upholstery** makes donor relaxed during blood donation.
- **Sleek design** of the donor couch reduce space constraints.
- **Interface** for blood collection monitor.

SPECIFICATIONS | SLEEK DONOR COUCH SC100

Power source	230 +/-10%VAC, 50/60Hz
Classification	Protection against electrical shock class I
Lifting capacity.	150 Kg(Max)
Weight	< 60 Kg
Movement / Position	Any Intermediate Position between sitting / Vasovagal attack position
Arm rest movement	Both swivellable and height adjustable.
Ingress of water	IPX0
Number of actuator	1
Control / Initiation of motor	Rocker switch
Power consumption	Operating : 120W, Standby:10W(max)
Dimensions	Length of backrest : 725 ±20 mm (2.4 ft) Length of seat : 1050 ±20 mm (3.5 ft) Seat width : 500 ±20 mm (1.7 ft) Seat height : 600 ±20 mm (2.0 ft) Armrest length : 500 ±20 mm (1.7 ft) Armrest width : 150 ±25 mm (0.5 ft) Overall length (trendelenburg position) : 1700 ±20mm Overall width (both armrests fully open) : 1530 ±20mm Overall height (sitting position) : 1200 ±20mm
Upholstery	Soft upholstery of 4 inch thickness
Maneuverability	Provided with lockable castor wheels/clamps for mounting in van
Accessories	IV Stand, Paper roll holder, BCM interface
Manufacturing standard	ISO 9001: 2008 , EN ISO 13485 : 2012
Standards & approvals	CE mark, EUroHS
Medical Device Classification	Class I

REGISTERED OFFICE

TERUMO PENPOL Private Limited
I-2, JAWAHAR NAGAR
KOWDIAR P.O, TRIVANDRUM – 695003
KERALA,INDIA.
Ph.No : 0471-3015500 / 501.
Fax: 0471-2721519
CIN: U33112KL1985PTC004531

www.terumopenpol.com
www.terumobct.com
pmt@terumopenpol.com

Toll Free No : 1800 3000 7070

OUR PRODUCT RANGE

BLOOD BAGS

DONOR STATION

BLOOD COLLECTION MONITOR

TUBE SEALER

PORTABLE TUBE SEALER

DUAL HEAD TUBE SEALER

COMPOSCALE + INTELLIGENT EXPRESSOR

CRYO BATH

BLOOD BANK REFRIGERATOR

PLATELET AGITATOR AND INCUBATOR

DEEP FREEZERS

PLASMA BATH

FILTERS*

BLOOD TRANSFUSION SET*

AUTOMATIC COMPONENT EXTRACTOR*

STERILE CONNECTING DEVICE*

APHERESIS SYSTEM AND KITS*

* Additional Products for Indian market only.
Products from Terumo Corporation / Terumo BCT

TERUMO PENPOL®