

## **BLOOD AND FLUID WARMER**

### **MODEL: BW 20 MBPL**

**-ALARM INDICATOR FOR LOW AND HIGH TEMPERATURE.**

**-DURABILITY AND EASY OF USE.**

**-USER-FRIENDLY**



**MEHROTRA BIOTECH PVT. LTD.**

M&M TOWER B-22 23 VIBHUTI KHAND GOMTI NAGAR,

Lucknow 226010

Tel + 91 522- 2720663

E-mail: sales@mehrotrabiotech.com, contact@mehrotrabiotech.com

Website: www.mehrotrabiotech.com

**ISO 9001: 2015**  
**ISO 13485: 2016**



A **Blood and Fluid** is Medical equipment used in the healthcare industry for warming blood products and Fluids at 30°C to 45 °C. Standard temperature for warming blood products is 37 °C and warming cabinets maintain the accuracy of typically +/- 0.1 °C.

Blood and Fluid warming cabinets are equipped with a fully automatic temperature controller and USB port. Neoprene Gasket saves the cabinet from any kind of leakage and temperature loss. Our warming cabinets are compact, lightweight & user-friendly.

#### TECHNICAL SPECIFICATIONS:

- External case treated with a corrosive resistant epoxy coating.
- Priming Volume in mL: <=50.
- Display: LCD
- Technology of heat exchanger: Dry Heat.
- Heating elements in the ribs (at the bottom) and three sides
- Temperature ranging from 30°C to 45 °C with 0.5 °C sensitivity.
- High Temp cut off in 41°C -43°C.
- Warm-up time of heat exchanger in min: <=3
- Input Power Supply: 230 V, 50Hz A.C., Single phase.
- CE, ISO 9001:2015, ISO 13485:2016 certified.

#### MEHROTRA BIOTECH PRIVATE LIMITED

**HEAD OFFICE: LUCKNOW ADDRESS** -M&M TOWER B-22 23 VIBHUTI KHAN GOMTI NAGAR, Lucknow 226010  
Tel + 91 522- 2720663 Fax + 91 522- 4004804.

**BRANCH OFFICE: DELHI ADDRESS**-302, Building No.38, 3rd floor Ansal Tower, Nehru Place New Delhi-110019  
Phone number- +91 9654785888, +91-011 41021495.

**BRANCH OFFICE: HYDERABAD ADDRESS**-G1, 6-3-1177, A/12, HUSSAIN NAGAR, KUNDANBAUG , HYDERABAD 500016.  
Phone number- 09911675755, 08447143644

**BRANCH OFFICE: MUMBAI ADDRESS**- 301, SUDAMA SHILPA, DR. R. P. ROAD DATTANAGAR, DOMBIVILLI (EAST)-421201,  
MUMBAI.  
Phone number: 9769173005

E-mail: [sales@mehrotrabiotech.com](mailto:sales@mehrotrabiotech.com), [service@mehrotrabiotech.com](mailto:service@mehrotrabiotech.com),

Website: [www.mehrotrabiotech.com](http://www.mehrotrabiotech.com).

ISO 9001: 2015  
ISO 13485: 2016

