



Front View



Back View

PORTABLE DONOR COUCH

- DC-101 is a portable Blood Collection Donor Chair which is easy to carry for any blood donation position. Move manually arm chair to full recliner with a comfortable blood donation condition.
- Chair can be fold and can be transported anywhere .
- Zero gravity reclining.
- The Poly Vinyl cushion make this chair Ultra Comfortable & Solid stainless Steel makes it extra strong & screwless riveted to reduce the risk of frequent mobility.

Specifications

- **Material used:** Stainless Steel AISI 304 ,16 gauge

Dimensions:-

- Length of back rest :34 +/-2 cm
- Width : 48 +/- 2cm
- Length of arm rest :46 +/-2 cm
- Width of arm rest :16 +/-2 cm
- Height of seat :64 +/- 2 cm

Mehrotra Biotech Pvt. Ltd.

M&M TOWER, B22-23, VIBHUTI KHAND, GOMTI NAGAR, LUCKNOW

Tel + 91 522- 4004768, 2326118 Fax + 91 522- 4004804

E-mail: devesh@mehrotrabiotech.com,

Website: www.mehrotrabiotech.com

ISO 9001: 2015
ISO 13485: 2016



TECHNICAL SPECIFICATION

- Vasovagal attack position height: 78 +/- 2 cm
- Max. Height of back rest: 120 +/- 2 cm
- Lifting Capacity:-155 Kg
- Weight:-15 kg with all accessories 20 kg.
- Arm Rest:-Both Left/right, and height adjustable
- Operation:-Manual

ACCESSORIES:

Comes with a Trolley which can be used for transportation of Donor Couch.



Two Donor Couches packed in a Trolley

Applications: Portable Donor Couch is important for blood donation camps which are organized by Blood Banks at some other/remote locations far from the blood center. This device is mandatory by Drug Controller of India and it also protects blood donor from any kind of mis-happening during Blood donation.

Mehrotra Biotech Pvt. Ltd.

M&M TOWER, B22-23, VIBHUTI KHAND, GOMTI NAGAR, LUCKNOW

Tel + 91 522- 4004768, 2326118 Fax + 91 522- 4004804

E-mail: devesh@mehrotrabiotech.com,

Website: www.mehrotrabiotech.com

ISO 9001: 2015
ISO 13485: 2016

