


INCUBATOR

- 1. PRECISION THERMOSTAT
CONTROLLED**
- 2. DURABLE CONSTRUCTION**
- 3. PRECISE TEMPERATURE
ACCURACY & UNIFORMITY**
- 4. DUAL DOOR**
- 5. ADJUSTABLE PERFORATED
SHELVES**



MEHROTRA BIOTECH PRIVATE LIMITED

M&M TOWER, B-22 23, VIBHUTI KHAND, GOMTI NAGAR,
Lucknow 226010.

Tel + 91 522- 2720663

E-mail: sales@mehrotrabiotech.com,

Website: www.mehrotrabiotech.com



TECHNICAL SPECIFICATION:

- Temperature ranging from Ambient + 5°C to 100°C.
- Temperature accuracy: ±1°C @37°C
- Material of construction (External): Corrosive resistant MS Powder coated.
- Material of construction (Internal): Stainless Steel AISI 304.
- Double walled construction.
- Stainless steel shelf (adjustable).
- Dual Door with Inner door made up of acrylic.
- All around seal on the doors using synthetic rubber gasket.
- Triple side heating.
- Adjustable air ventilators on both the sides.
- Precision thermostat controlled.
- Door open alarm present.
- Air circulation by natural convection.
- Power Supply: 230V, 50Hz.

STANDARD SIZES/ DIMENSIONS:


MODEL	DIMENSION (INCH)	VOLUME (approx)	SHELVES	WATT
MI 30	12x12x12	30 L	2	250 W
MI 50	14x14x14	50 L	2	250 W
MI 60	15x15x16	60 L	2	280 W
MI 90	18x18x18	90 L	2	400 W
MI 120	18x18x24	120 L	3	400 W
MI 150	21x21x21	150 L	3	500 W
MI 200	24x24x24	200 L	3	600 W

OPTIONAL:

- Fan Assisted Air Circulation.
- 21 CFR part 11 Software.
- 5” TFT Touch Screen Controller.

Display type:

-Digital Display



- Double Segment Digital Display.
- Microprocessor controlled.
- Audio-Visual alarm of Temperature with set, process value.

-TFT Touch Screen Display



- 5" TFT Touch Screen controller.
- Microprocessor controlled.
- Visual and Acoustic alarm for Door Open, Power Failure, High/Low Temperature, Low battery.
- Digital door micro-switch.
- Keypad Lock for no data manipulation.
- Data Logging.
- USB port.

MEHROTRA BIOTECH PRIVATE LIMITED

ADDRESS -M&M TOWER, B-22 23, Vibhuti Khand, Gomti Nagar, Lucknow-226010.

Tel: 0522- 2720663.

E-mail: sales@mehrotrabiotech.com, service@mehrotrabiotech.com,

Website: www.mehrotrabiotech.com.