

MEHROTRA GROUP



**HOUSE OF BIOMEDICAL RESEARCH
& SCIENTIFIC EQUIPMENT**

MBPE-II MANUAL PLASMA EXPRESSOR




➤APPLICATION

Plasma Expressor is used for separating plasma in blood banks, hospitals, and other facilities. After the blood inside the blood bag is layered, the pressure of the plasma extractor can push out the different components of the blood from top to bottom.

Model	MBPE-II
Applicable Blood Bag Capacity	100 ~ 400ml
Clamping Force	20N
Material	Stainless steel
Dimensions(W×D×H)(mm)	250x165x250
Shipping size(W×D×H)(mm)	320x320x270
N.W./G.W. (kg)	2.5/3.1

MEHROTRA GROUP



HOUSE OF BIOMEDICAL RESEARCH & SCIENTIFIC EQUIPMENT

M&M TOWER B-22 23, VIBHUTI KHAND, GOMTI NAGAR, LUCKNOW- 226010

Tel + 91 522- 2720663, 9838203524,9838203399,9839014323

E-mail: sales@mehrotrabiotech.com, service@mehrotrabiotech.com.

Website: www.mehrotrabiotech.com

OUR BRANCH OFFICES

LUCKNOW

DELHI

MUMBAI

HYDERABAD

BANGALORE

KOLKATA

AHMEDABAD

