

MEHROTRA BIOTECH-HOT PLATE

FEATURES

- Unique design made of perforated PCRC heavy duty sheet duly epoxy powder coated with gap of cast iron top from the body.
- Kanthal A-1 superior quality coil shaped heater duly embedded in the cast iron plate to avoid any heat loss having Cera Blanket insulation.
- All electrical connections are fire proof having sleeve for Heat/Temperature Protection.
- Temp. Range is ambient to 3500° C.
- Temperature controlled by SUNVIC Energy Regulator having marked 0-100, at minimum position it release minimum energy and at max. Position it release maximum energy.

OPTIONAL:

- Thermostat control.
- Pyrometer with thermocouple
Solid State Digital Temperature indicator-cum-controller.



MEHROTRA BIOTECH PVT. LTD.

M&M TOWER, B-22, 23,
VIBHUTI KHAND, GOMTI NAGAR,
LUCKNOW.

Tel + 91 522- 2720663.

E-mail: sales@mehrotrabiotech.com, contact@mehrotrabiotech.com,

Website: www.mehrotrabiotech.com.